

Medela Clario Toni suction pump

Pump Specifications																																	
Item numbers	<p>AC/DC version</p> <table border="0"> <tr> <td>014.0022</td> <td>100–240V/50–60Hz</td> <td>AUS/NZ plug</td> <td>A</td> </tr> <tr> <td>014.0023</td> <td>100–240V/50–60Hz</td> <td>Euro plug</td> <td>A</td> </tr> <tr> <td>014.0024</td> <td>100–240V/50–60Hz</td> <td>Euro plug</td> <td>B</td> </tr> <tr> <td>014.0025</td> <td>100–240V/50–60Hz</td> <td>UK Plug</td> <td>A</td> </tr> </table> <p>Deluxe version (AC/DC + SoftCase)</p> <table border="0"> <tr> <td>014.0026</td> <td>100–240V/50–60Hz</td> <td>AUS/NZ plug</td> <td>A</td> </tr> <tr> <td>014.0027</td> <td>100–240V/50–60Hz</td> <td>Euro plug</td> <td>A</td> </tr> <tr> <td>014.0028</td> <td>100–240V/50–60Hz</td> <td>Euro plug</td> <td>B</td> </tr> <tr> <td>014.0029</td> <td>100–240V/50–60Hz</td> <td>UK Plug</td> <td>A</td> </tr> </table>	014.0022	100–240V/50–60Hz	AUS/NZ plug	A	014.0023	100–240V/50–60Hz	Euro plug	A	014.0024	100–240V/50–60Hz	Euro plug	B	014.0025	100–240V/50–60Hz	UK Plug	A	014.0026	100–240V/50–60Hz	AUS/NZ plug	A	014.0027	100–240V/50–60Hz	Euro plug	A	014.0028	100–240V/50–60Hz	Euro plug	B	014.0029	100–240V/50–60Hz	UK Plug	A
014.0022	100–240V/50–60Hz	AUS/NZ plug	A																														
014.0023	100–240V/50–60Hz	Euro plug	A																														
014.0024	100–240V/50–60Hz	Euro plug	B																														
014.0025	100–240V/50–60Hz	UK Plug	A																														
014.0026	100–240V/50–60Hz	AUS/NZ plug	A																														
014.0027	100–240V/50–60Hz	Euro plug	A																														
014.0028	100–240V/50–60Hz	Euro plug	B																														
014.0029	100–240V/50–60Hz	UK Plug	A																														
Instructions for use	<table border="0"> <tr> <td>177.1405</td> <td>EN, DE, FR, IT, NL, SV, FI, NO, ES, PT (A)</td> </tr> <tr> <td>177.1406</td> <td>ET, HU, CS, LT, LV, PL, SL, SK, DA, EL (B)</td> </tr> </table>	177.1405	EN, DE, FR, IT, NL, SV, FI, NO, ES, PT (A)	177.1406	ET, HU, CS, LT, LV, PL, SL, SK, DA, EL (B)																												
177.1405	EN, DE, FR, IT, NL, SV, FI, NO, ES, PT (A)																																
177.1406	ET, HU, CS, LT, LV, PL, SL, SK, DA, EL (B)																																
Service manual	200.4001 EN, DE, FR																																
<p>Max. vacuum [kPa/mmHg] Measured at 0 m, atmospheric pressure: 1013.25 hPa Please note: vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).</p>	<p>Tolerance:</p> <table border="0"> <tr> <td>min. –8/–60</td> <td>± 15%</td> </tr> <tr> <td>med. –13/–100</td> <td>± 15%</td> </tr> <tr> <td>max. –80/–600</td> <td>– 15%</td> </tr> </table>	min. –8/–60	± 15%	med. –13/–100	± 15%	max. –80/–600	– 15%																										
min. –8/–60	± 15%																																
med. –13/–100	± 15%																																
max. –80/–600	– 15%																																
Flow [l/min]	15 l/min. ± 10%																																
Dimensions [mm]	223 x 255 x 95																																
Weight [kg]	2.0																																
Drive unit	Membrane																																
Power supply	100–240V, 50–60 Hz																																
Power consumption [V]	12V, 3.4 A																																
Noise level [dB(A)]; measured in Medela acoustic room at 1 m distance	64.3																																
Protection against electric shock	Type BF																																
Medical Device Classification – According to 93/42/EEC	Class IIa																																
Protection against electric shock	Class II Equipment																																
Ingress Protection	IP21																																
Battery	NiMH																																
Operation from battery (at full charge)	>50 minutes																																
Charging time of battery	5 hours																																
Overflow protection	Threefold hygienic design: – Floater in canister prevents overflow – SafetyChamber – Membrane system No overflow can occur into the pump																																
Warranty [years]	2; Battery 6 months																																
Operating temperature, relative humidity, pressure	5 to 40 °C, 15–93 %, 70–106 kPa																																
Transport/storage temperature, relative humidity, pressure	–20 to +50 °C, 15–93 %, 70–106 kPa																																
Made in	Switzerland																																

Accessories	Item number
Canister set	074.0001
SafetyChamber Clario Toni	074.0030
Disposable Fingertip tubing	074.0003
Car connection cable	074.0017
SoftCase carrying bag	074.0007
Vacuum gauge	074.0006

CAUTION: Federal US law restricts this device to sale by or on the order of a physician.

Medical Vacuum Technology
for Healthcare Professionals

**Please contact us or your local Medela
representative for details.**

 Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

Local contact: