

Tubing Portfolio

SAFE AND RELIABLE CONNECTION TO MEDELA DEVICES



Tested compatibility and safety

- Seamless compatibility and efficient functionality with every Medela suction device

Many choices





- Extensive range of Medela tubing solutions
- Large choice of material and lengths

Many proven choices

Annex I of Directive 93/42/EEC for medical devices states: "If the device is intended for use in combination with other devices or equipment, the whole combination, including the connection system, must be safe and must not impair the specified performances of the devices." Every component in the Medela suction family (pumps, collection systems, and tubing) has been tested in combination with other Medela suction components, and has been proven safe, reliable and fully compliant with the Directive, as part of a Medela-compatible system. Medela offers a variety of components and configurations to choose from, according to your needs:


Vacuum tubing (from suction pump to collection system)

Vacuum tubing, silicone, ID x OD: 7x12 mm



	077.0054: 50cm 077.0902: 200cm 077.0900: 2500cm	077.0901: 100cm 077.0905: 500cm
	077.0931: 25 cm 077.0185: 60 cm	(1 st choice for Vario)
	077.0921: 30 cm 077.0922: 60 cm 077.0189: 300 cm	(1 st choice for Basic/Dominant Flex)
	077.0154: 25 cm 077.0912: 60 cm	077.0911: 33 cm 077.0913: 100 cm

Patient tubing (from collection system to patient applied part)



Disposable, ID x OD: 7x10 mm

	077.0170: 150 cm, sterile
---	---------------------------

Disposable, ID x OD: 6x8 mm

	077.0184: 150 cm, sterile
	077.0965: 150 cm, sterile


Disposable, ID x OD: 6x8 mm

	077.0951: 180 cm, sterile
	077.0952: 180 cm



Disposable, ID x OD: 12x18 mm

	077.0963: 180 cm for large lid 077.0430 or direct connection on liner
---	---





Silicone, ID x OD: 5x10 mm

	077.0970: 2500 cm
---	-------------------


Silicone, ID x OD: 7x12 mm


	077.0053: 32 cm 077.0056: 150 cm 077.0164: 300 cm	077.0055: 120 cm 077.0051: 200 cm 077.0052: 2500 cm
	077.0961: 180 cm, silicone (green)	

Tubing for high volume setups


	077.0563 Change-over valve with vacuum tubing to manually switch vacuum from one jar to the other
	077.0095: 25 cm, serial connection, vacuum tubing, DCS
	077.0158: 34 cm (100 pc) serial connection DCS, disposable ID x OD: 8 x 12 mm 077.0160: 34 cm serial connection DCS, disposable ID x OD: 8 x 12 mm
	077.0542: 50 cm serial connection RCS, silicone ID x OD: 7 x 12 mm

Tested compatibility with:

	Vario		
	18	18c/i	8*
Max. vac. kPa	-75	-55	-9
Flow l/min.	18	18	8
* Vario 8 to be used with Reusable Collection System (RCS) only			

	Dominant Flex		Basic
Max. vac. kPa	-95		-90
Flow l/min.	40, 50, 60		30

	Medela Disposable Collection System (DCS) Sizes: 1.5L / 2.5L with/without solidifier
---	---

	Medela Reusable Collection System (RCS) Sizes: 0.25, 0.5, 1, 2, 3, 5L
---	---

Medical Vacuum Technology for
Healthcare Professionals

Please contact us or your local
Medela Representative for details.

 Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

Local contact

CE 0123